

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Eotaxin Antibody |
| Catalog Number: | C30181 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from the Internal region of human CCL11. AA range:41-90 |
| Concentration: | 1 mg/ml |
| Gene Name: | CCL11 |
| Protein Name: | Eotaxin |
| Other Name: | CCL11; SCYA11; Eotaxin; C-C motif chemokine 11; Eosinophil chemotactic protein; Small-inducible cytokine A11 |
| Human Swiss Prot No.: | P51671 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Eotaxin Polyclonal Antibody detects endogenous levels of Eotaxin protein. |
| Storage and Stability: | -20°C/1 year |