

Antibodies Product Datasheet

Product Name:	COX4 Antibody
Catalog Number:	C47463
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:1000 IHC: 1:500-1000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of COX4
Concentration:	
Gene Name:	COX4I1
Protein Name:	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial
Other Name:	COX4I1; COX4; Cytochrome c oxidase subunit 4 isoform 1, mitochondrial; Cytochrome c oxidase polypeptide IV; Cytochrome c oxidase subunit IV isoform 1; COX IV-1; YM3351
Human Swiss Prot No.:	P13073
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous COX4 protein.
Storage and Stability:	-20°C/1 year