

Antibodies Product Datasheet

Product Name:	FA11 (heavy chain, Cleaved-Glu19) Antibody
Catalog Number:	L0423
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FA11 (heavy chain, Cleaved-Glu19)
Concentration:	1 mg/ml
Gene Name:	F11
Protein Name:	FA11 (heavy chain, Cleaved-Glu19)
Other Name:	Coagulation factor XI (FXI) (EC 3.4.21.27) (Plasma thromboplastin antecedent) (PTA) [Cleaved into: Coagulation factor XIa heavy chain; Coagulation factor XIa light chain]
Human Swiss Prot No.:	P03951
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse FA11 (heavy chain, Cleaved-Glu19, protein was cleaved amino acid sequence between 18-19)
Storage and Stability:	-20°C/1 year