

Antibodies Product Datasheet

Product Name:	p23 Antibody
Catalog Number:	B0583
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TEBP. AA range:79-128
Concentration:	1 mg/ml
Gene Name:	PTGES3
Protein Name:	Prostaglandin E synthase 3
Other Name:	PTGES3; P23; TEBP; Prostaglandin E synthase 3; Cytosolic prostaglandin E2 synthase; cPGES; Hsp90 co-chaperone; Progesterone receptor complex p23; Telomerase-binding protein p23
Human Swiss Prot No.:	Q15185
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	p23 Polyclonal Antibody detects endogenous levels of p23 protein.
Storage and Stability:	-20°C/1 year