

Antibodies Product Datasheet

Product Name:	GABAB R1 Antibody
Catalog Number:	G266
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GABBR1. AA range:891-940
Concentration:	1 mg/ml
Gene Name:	GABBR1
Protein Name:	Gamma-aminobutyric acid type B receptor subunit 1
Other Name:	GABBR1; GPRC3A; Gamma-aminobutyric acid type B receptor subunit 1; GABA-B receptor 1; GABA-B-R1; GABA-BR1; GABABR1; Gb1
Human Swiss Prot No.:	Q9UBS5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABAB R1 Polyclonal Antibody detects endogenous levels of GABAB R1 protein.
Storage and Stability:	-20°C/1 year