

Antibodies Product Datasheet

Product Name:	Presenilin 1 (Phospho-Ser357) Antibody
Catalog Number:	A8404
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PSEN1 around the phosphorylation site of Ser357. AA range:323-372
Concentration:	1 mg/ml
Gene Name:	PSEN1
Protein Name:	Presenilin-1
Other Name:	PSEN1; AD3; PS1; PSNL1; Presenilin-1; PS-1; Protein S182
Human Swiss Prot No.:	P49768
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Presenilin 1 (S357) Polyclonal Antibody detects endogenous levels of Presenilin 1 protein only when phosphorylated at S357.
Storage and Stability:	-20°C/1 year