

Antibodies Product Datasheet

Product Name:	Cleaved-C1r HC (R463) Antibody
Catalog Number:	L0229
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human C1R. AA range:414-463
Concentration:	1 mg/ml
Gene Name:	C1R
Protein Name:	Complement C1r subcomponent
Other Name:	C1R; Complement C1r subcomponent; Complement component 1 subcomponent r
Human Swiss Prot No.:	P00736
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-C1r HC (R463) Polyclonal Antibody detects endogenous levels of fragment of activated C1r HC protein resulting from cleavage adjacent to R463.
Storage and Stability:	-20°C/1 year