

## Antibodies Product Datasheet

<b>Product Name:</b>	TIP60 Monoclonal Antibody
<b>Catalog Number:</b>	N0619
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,IHC-p,IF(paraffin section),ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of human TIP60 expressed in E. Coli.
<b>Concentration:</b>	
<b>Gene Name:</b>	KAT5
<b>Protein Name:</b>	Histone acetyltransferase KAT5
<b>Other Name:</b>	KAT5; HTATIP; TIP60; Histone acetyltransferase KAT5; 60 kDa Tat- interactive protein; Tip60; Histone acetyltransferase HTATIP; HIV-1 Tat interactive protein; Lysine acetyltransferase 5; cPLA(2)-interacting protein
<b>Human Swiss Prot No.:</b>	Q92993
<b>Formulation:</b>	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
<b>Specificity:</b>	TIP60 Monoclonal Antibody detects endogenous levels of TIP60 protein.
<b>Storage and Stability:</b>	-20°C/1 year