

Antibodies Product Datasheet

Product Name:	PGPEP1 Antibody
Catalog Number:	C50346
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PGPEP1
Concentration:	1 mg/ml
Gene Name:	PGPEP1 PGPI
Protein Name:	Pyroglutamyl-peptidase 1 (EC 3.4.19.3) (5-oxopropyl-peptidase) (Pyroglutamyl aminopeptidase I) (PAP-I) (Pyroglutamyl-peptidase I) (PGP-I) (Pyrrolidone-carboxylate peptidase)
Other Name:	
Human Swiss Prot No.:	Q9NXJ5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PGPEP1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year