

Antibodies Product Datasheet

Product Name:	PLA1A Antibody
Catalog Number:	C17709
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLA1A. AA range:381-430
Concentration:	1 mg/ml
Gene Name:	PLA1A
Protein Name:	Phospholipase A1 member A
Other Name:	PLA1A; NMD; PSPLA1; Phospholipase A1 member A; Phosphatidylserine-specific phospholipase A1; PS-PLA1
Human Swiss Prot No.:	Q53H76
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PLA1A Polyclonal Antibody detects endogenous levels of PLA1A protein.
Storage and Stability:	-20°C/1 year