

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Bcl-2 Monoclonal Antibody(6B5)
<b>Catalog Number:</b>	N1305
<b>Reactivity:</b>	Human,Mouse,Rat,chicken
<b>Applications:</b>	IF/ICC,WB,IHC-p
<b>Dilutions:</b>	IF: 1:50-200 WB: 1:1000~2000 IHC: 1:200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of Bcl-2
<b>Concentration:</b>	
<b>Gene Name:</b>	BCL2
<b>Protein Name:</b>	Apoptosis regulator Bcl-2
<b>Other Name:</b>	BCL2; Apoptosis regulator Bcl-2
<b>Human Swiss Prot No.:</b>	P10415
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous Bcl-2 proteins.
<b>Storage and Stability:</b>	-20°C/1 year