
Antibodies Product Datasheet

Product Name:	Brn-3 Antibody
Catalog Number:	C13491
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human POU4F3. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	POU4F3
Protein Name:	POU domain class 4 transcription factor 3
Other Name:	POU4F3; BRN3C; POU domain; class 4, transcription factor 3; Brain-specific homeobox/POU domain protein 3C; Brain-3C; Brn-3C
Human Swiss Prot No.:	Q15319
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Brn-3 Polyclonal Antibody detects endogenous levels of Brn-3 protein.
Storage and Stability:	-20°C/1 year