

Antibodies Product Datasheet

Product Name:	DGK- ϵ Antibody
Catalog Number:	C10697
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DGKE. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	DGKE
Protein Name:	Diacylglycerol kinase epsilon
Other Name:	DGKE; DAGK5; Diacylglycerol kinase epsilon; DAG kinase epsilon; Diglyceride kinase epsilon; DGK-epsilon
Human Swiss Prot No.:	P52429
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DGK- ϵ Polyclonal Antibody detects endogenous levels of DGK- ϵ protein.
Storage and Stability:	-20°C/1 year