

Antibodies Product Datasheet

Product Name:	Integrin β 3 Antibody
Catalog Number:	C30370
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human ITGB3. AA range:621-670
Concentration:	1 mg/ml
Gene Name:	ITGB3
Protein Name:	Integrin beta-3
Other Name:	ITGB3; GP3A; Integrin beta-3; Platelet membrane glycoprotein IIIa; GPIIIa; CD61
Human Swiss Prot No.:	P05106
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Integrin β 3 Polyclonal Antibody detects endogenous levels of Integrin β 3 protein.
Storage and Stability:	-20°C/1 year