

Antibodies Product Datasheet

Product Name:	Rag A/B Antibody
Catalog Number:	C18211
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RRAGA/B. AA range:264-313
Concentration:	1 mg/ml
Gene Name:	RRAGA/RRAGB
Protein Name:	Ras-related GTP-binding protein A/B
Other Name:	RRAGA; Ras-related GTP-binding protein A; Rag A; RagA; Adenovirus E3 14.7 kDa-interacting protein 1; FIP-1; RRAGB; Ras-related GTP-binding protein B; Rag B; RagB
Human Swiss Prot No.:	Q7L523/Q5VZM2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rag A/B Polyclonal Antibody detects endogenous levels of Rag A/B protein.
Storage and Stability:	-20°C/1 year