
Antibodies Product Datasheet

Product Name:	PKC δ Antibody
Catalog Number:	B7196
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PKC delta. AA range:612-661
Concentration:	1 mg/ml
Gene Name:	PRKCD
Protein Name:	Protein kinase C delta type
Other Name:	PRKCD; Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta
Human Swiss Prot No.:	Q05655
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PKC δ Polyclonal Antibody detects endogenous levels of PKC δ protein.
Storage and Stability:	-20°C/1 year