

Antibodies Product Datasheet

Product Name:	PKC δ (Phospho-Thr505) Antibody
Catalog Number:	A9464
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human PKC δ (Thr505)
Concentration:	1 mg/ml
Gene Name:	PRKCD
Protein Name:	PKC δ (Thr505)
Other Name:	Protein kinase C delta type (EC 2.7.11.13) (Tyrosine-protein kinase PRKCD) (EC 2.7.10.2) (nPKC-delta)
Human Swiss Prot No.:	Q05655
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat PKC δ (phospho-Thr505)
Storage and Stability:	-20°C/1 year