
Antibodies Product Datasheet

Product Name:	ATG13 Antibody
Catalog Number:	C42080
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ATG13
Concentration:	1 mg/ml
Gene Name:	ATG13 KIAA0652
Protein Name:	ATG13
Other Name:	Autophagy-related protein 13; YM33059
Human Swiss Prot No.:	O75143
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.19% sodium azide.
Specificity:	This antibody detects endogenous levels of ATG13 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year