

Antibodies Product Datasheet

Product Name:	SCAF1 Antibody
Catalog Number:	C18936
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SFRS19. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	SCAF1
Protein Name:	Splicing factor arginine/serine-rich 19
Other Name:	SCAF1; SFRS19; SRA1; Splicing factor; arginine/serine-rich 19; SR-related and CTD-associated factor 1; SR-related-CTD-associated factor; SCAF; Serine arginine-rich pre-mRNA splicing factor SR-A1; SR-A1
Human Swiss Prot No.:	Q9H7N4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SCAF1 Polyclonal Antibody detects endogenous levels of SCAF1 protein.
Storage and Stability:	-20°C/1 year