## **Assay Biotechnology Company**

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

## **Antibodies Product Datasheet**

**Product Name:** Topo IIα (Phospho-Ser1525) Antibody

Catalog Number: A8439

Reactivity: Human, Mouse, Rat

**Applications:** WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000. Not yet tested in other applications.

Modification: Phospho Specific

**Source:** Polyclonal Rabbit

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Immunogen:** The antiserum was produced against synthesized peptide derived from

human TOP2A around the phosphorylation site of Ser1525. AA

range:1482-1531

Concentration: 1 mg/ml

Gene Name: TOP2A

**Protein Name:** DNA topoisomerase 2-alpha

Other Name: TOP2A; TOP2; DNA topoisomerase 2-alpha; DNA topoisomerase II;

alpha isozyme

**Human Swiss Prot No.:** P11388

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

**Specificity:** Phospho-Topo IIα (S1525) Polyclonal Antibody detects endogenous

levels of Topo IIα protein only when phosphorylated at S1525.

Storage and Stability: -20°C/1 year