

Antibodies Product Datasheet

Product Name:	Doublecortin (Phospho-Ser378) Antibody
Catalog Number:	A0637
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Doublecortin (Phospho-Ser378)
Concentration:	1 mg/ml
Gene Name:	DCX DBCN LISX
Protein Name:	Doublecortin (Phospho-Ser378)
Other Name:	Neuronal migration protein doublecortin (Doublin) (Lissencephalin-X) (Lis-X)
Human Swiss Prot No.:	O43602
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of Doublecortin (Phospho-Ser378) at Human:S378, Mouse:S297, Rat:S297
Storage and Stability:	-20°C/1 year