

Antibodies Product Datasheet

Product Name:	PAR4 Antibody
Catalog Number:	C0309
Reactivity:	Human, Mouse, Rat, Cow
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Prostate Apoptosis Response protein-4. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	PAWR
Protein Name:	PRKC apoptosis WT1 regulator protein
Other Name:	PAWR; PAR4; PRKC apoptosis WT1 regulator protein; Prostate
Human Swiss Prot No.:	apoptosis response 4 protein; Par-4 Q96IZ0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PAR4 Polyclonal Antibody detects endogenous levels of PAR4 protein.
Storage and Stability:	-20°C/1 year