

Antibodies Product Datasheet

Product Name:	PRMT6 mouse mAb
Catalog Number:	N0788
Reactivity:	Human
Applications:	WB,IHC,ICC
Dilutions:	wb 1:1000 icc 1:300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human PRMT6 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	prmt6
Protein Name:	
Other Name:	ANM6_HUMAN;Chromobox protein homolog 7;FLJ10559;FLJ51477;Heterogeneous nuclear ribonucleoprotein methyltransferase like protein 6;Heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 6;Histone-arginine N-methyltransferase PRMT6;HMT1 hnRN
Human Swiss Prot No.:	Q96LA8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PRMT6 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year