

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Syntaxin 1 (Phospho-Ser14) Antibody |
| Catalog Number: | A1238 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-p, IF/ICC, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Phospho Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human Syntaxin 1A around the phosphorylation site of Ser14. AA range:1-50 |
| Concentration: | 1 mg/ml |
| Gene Name: | STX1A |
| Protein Name: | Syntaxin-1A |
| Other Name: | STX1A; STX1; Syntaxin-1A; Neuron-specific antigen HPC-1 |
| Human Swiss Prot No.: | Q16623 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Phospho-Syntaxin 1 (S14) Polyclonal Antibody detects endogenous levels of Syntaxin 1 protein only when phosphorylated at S14. |
| Storage and Stability: | -20°C/1 year |