

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | Rabbit Anti Guinea Pig IgG (H&L) Biotin |
| Catalog Number: | SA0368 |
| Reactivity: | Guinea Pig |
| Conjugate: | Biotin |
| Applications: | ELISA, WB, IHC, ICC |
| Dilutions: | ELISA (1/5000 - 1/20000), WB (1/5000 - 1/20000), IHC (1/500 - 1/1000), ICC (1/500 - 1/1000) |
| Modification: | Not Modification Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Concentration: | 1mg/mL |
| Other Name: | Rabbit Anti-Guinea Pig IgG (H&L) Biotin |
| Formulation: | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide |
| Storage and Stability: | Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing |