

Antibodies Product Datasheet

Product Name:	KCNN4 (SK4) Antibody
Catalog Number:	C40062
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:50-100
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of KCNN4 (SK4)
Concentration:	1 mg/ml
Gene Name:	KCNN4
Protein Name:	Intermediate conductance calcium-activated potassium channel protein 4 (SK4) (SKCa 4) (SKCa4) (IKCa1) (IK1) (KCa3.1) (KCa4) (Putative Gardos channel)
Other Name:	KCNN4; IK1; IKCA1; KCA4; SK4; Intermediate conductance calcium-activated potassium channel protein 4; SK4; SKCa 4; SKCa4; IKCa1; IK1; KCa3.1; KCa4; Putative Gardos channel; YM3567
Human Swiss Prot No.:	O15554
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KCNN4(SK4) protein(A247) detects endogenous levels of KCNN4(SK4)
Storage and Stability:	-20°C/1 year