

Antibodies Product Datasheet

Product Name:	Caspase-10 B/C Antibody
Catalog Number:	C0135
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Caspase 10. AA range:430-479
Concentration:	1 mg/ml
Gene Name:	CASP10
Protein Name:	Caspase10
Other Name:	ALPS2; apoptotic protease Mch-4; CASP-10; CASPA; caspase-10; FADD-like ICE2; FAS-associated death domain protein interleukin-1B-convert
Human Swiss Prot No.:	Q92851-2/4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Caspase-10 B/C Polyclonal Antibody detects endogenous levels of Caspase-10 B/C protein.
Storage and Stability:	-20°C/1 year