

Antibodies Product Datasheet

Product Name:	CAR18 Antibody
Catalog Number:	C42868
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CAR18
Concentration:	1 mg/ml
Gene Name:	CARD18 ICEBERG UNQ5804/PRO19611
Protein Name:	CAR18
Other Name:	
Human Swiss Prot No.:	P57730
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CAR18 at Human
Storage and Stability:	-20°C/1 year