

Antibodies Product Datasheet

Product Name:	HSC70 Rabbit Antibody
Catalog Number:	C40081
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:500-2000,IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)
Concentration:	1 mg/ml
Gene Name:	HSPA8 HSC70 HSP73 HSPA10
Protein Name:	Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)
Other Name:	Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8); YM3715
Human Swiss Prot No.:	P11142
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous HSC70 protein
Storage and Stability:	-20°C/1 year