

---

## Antibodies Product Datasheet

<b>Product Name:</b>	CU082 Antibody
<b>Catalog Number:</b>	C40552
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	IHC-p, IF(paraffin section)
<b>Dilutions:</b>	IHC-p 1: 50-200
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human CU082
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	LINC00310 C21orf82 NCRNA00310
<b>Protein Name:</b>	CU082
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	P59036
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of CU082 at Human
<b>Storage and Stability:</b>	-20°C/1 year