

## Antibodies Product Datasheet

<b>Product Name:</b>	TOPBP1 (Phospho-Ser1159) Antibody
<b>Catalog Number:</b>	A9126
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human TOPBP1 (Phospho-Ser1159)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	TOPBP1 KIAA0259
<b>Protein Name:</b>	TOPBP1 (Phospho-Ser1159)
<b>Other Name:</b>	DNA topoisomerase 2-binding protein 1 (DNA topoisomerase II-beta-binding protein 1) (TopBP1) (DNA topoisomerase II-binding protein 1)
<b>Human Swiss Prot No.:</b>	Q92547
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.168% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of TOPBP1 (Phospho-Ser1159) at Human, Mouse, Rat
<b>Storage and Stability:</b>	-20°C/1 year