

Antibodies Product Datasheet

Product Name:	PKC ϵ Antibody
Catalog Number:	B0802
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PKC epsilon. AA range:688-737
Concentration:	1 mg/ml
Gene Name:	PRKCE
Protein Name:	Protein kinase C epsilon type
Other Name:	PRKCE; PKCE; Protein kinase C epsilon type; nPKC-epsilon
Human Swiss Prot No.:	Q02156
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PKC ϵ Polyclonal Antibody detects endogenous levels of PKC ϵ protein.
Storage and Stability:	-20°C/1 year