
Antibodies Product Datasheet

Product Name:	Dorfin Antibody
Catalog Number:	C15578
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RNF19A. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	RNF19A
Protein Name:	E3 ubiquitin-protein ligase RNF19A
Other Name:	RNF19A; RNF19; E3 ubiquitin-protein ligase RNF19A; Double ring-finger protein; Dorfin; RING finger protein 19A; p38
Human Swiss Prot No.:	Q9NV58
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dorfin Polyclonal Antibody detects endogenous levels of Dorfin protein.
Storage and Stability:	-20°C/1 year