

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | EphB4 Monoclonal Antibody |
| Catalog Number: | N0233 |
| Reactivity: | Human |
| Applications: | WB,IF/ICC,ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| Purification: | Affinity purification |
| Immunogen: | Purified recombinant extracellular fragment of human EphB4 fused with hlgGfc tag expressed in HEK293 cell line. |
| Concentration: | |
| Gene Name: | EPHB4 |
| Protein Name: | Ephrin type-B receptor 4 |
| Other Name: | EPHB4; HTK; MYK1; TYRO11; Ephrin type-B receptor 4; Hepatoma transmembrane kinase; Tyrosine-protein kinase TYRO11 |
| Human Swiss Prot No.: | P54760 |
| Formulation: | Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol. |
| Specificity: | EphB4 Monoclonal Antibody detects endogenous levels of EphB4 protein. |
| Storage and Stability: | -20°C/1 year |