

Antibodies Product Datasheet

Product Name:	ACBD3 Antibody
Catalog Number:	C50901
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ACBD3
Concentration:	1 mg/ml
Gene Name:	ACBD3 GCP60 GOCAP1 GOLPH1
Protein Name:	Golgi resident protein GCP60 (Acyl-CoA-binding domain-containing protein 3) (Golgi complex-associated protein 1) (GOCAP1) (Golgi phosphoprotein 1) (GOLPH1) (PBR- and PKA-associated protein 7) (Peripheral benzodiazepine receptor-associated protein PAP7)
Other Name:	
Human Swiss Prot No.:	Q9H3P7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ACBD3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year