

Antibodies Product Datasheet

Product Name:	DUSP5 Antibody
Catalog Number:	C51065
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DUSP5
Concentration:	1 mg/ml
Gene Name:	DUSP5 VH3
Protein Name:	Dual specificity protein phosphatase 5 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH3)
Other Name:	
Human Swiss Prot No.:	Q16690
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DUSP5 at Human,Rat
Storage and Stability:	Store at -20° C/1 year