

Antibodies Product Datasheet

Product Name:	ABHD6 Antibody
Catalog Number:	C14218
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABHD6. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	ABHD6
Protein Name:	Monoacylglycerol lipase ABHD6
Other Name:	ABHD6; Monoacylglycerol lipase ABHD6; 2-arachidonoylglycerol hydrolase; Abhydrolase domain-containing protein 6
Human Swiss Prot No.:	Q9BV23
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABHD6 Polyclonal Antibody detects endogenous levels of ABHD6 protein.
Storage and Stability:	-20°C/1 year