

Antibodies Product Datasheet

Product Name:	Wee1 (Phospho-Ser642) Antibody
Catalog Number:	A1186
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human WEE1 around the phosphorylation site of Ser642. AA range:597-646
Concentration:	1 mg/ml
Gene Name:	WEE1
Protein Name:	Wee1-like protein kinase
Other Name:	WEE1; Wee1-like protein kinase; WEE1hu; Wee1A kinase
Human Swiss Prot No.:	P30291
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Wee1 (S642) Polyclonal Antibody detects endogenous levels of Wee1 protein only when phosphorylated at S642.
Storage and Stability:	-20°C/1 year