

## Antibodies Product Datasheet

<b>Product Name:</b>	IGF2R (Phospho-Ser2484) Antibody
<b>Catalog Number:</b>	A1051
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC-p, IF(paraffin section), WB
<b>Dilutions:</b>	IHC-p 1:50-200, WB 1:500-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human IGF2R (Phospho-Ser2484)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	IGF2R MPRI
<b>Protein Name:</b>	IGF2R (Phospho-Ser2484)
<b>Other Name:</b>	Cation-independent mannose-6-phosphate receptor (CI Man-6-P receptor) (CI-MPR) (M6PR) (300 kDa mannose 6-phosphate receptor) (MPR 300) (Insulin-like growth factor 2 receptor) (Insulin-like growth factor II receptor) (IGF-II receptor) (M6P/IGF2 receptor) (
<b>Human Swiss Prot No.:</b>	P11717
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous phospho levels of IGF2R (Phospho-Ser2484) at Human:S2484, Mouse:S2476, Rat:S2472
<b>Storage and Stability:</b>	-20°C/1 year