

Antibodies Product Datasheet

Product Name:	BRS-3 Antibody
Catalog Number:	G034
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BRS3. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	BRS3
Protein Name:	Bombesin receptor subtype-3
Other Name:	BRS3; Bombesin receptor subtype-3; BRS-3
Human Swiss Prot No.:	P32247
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BRS-3 Polyclonal Antibody detects endogenous levels of BRS-3 protein.
Storage and Stability:	-20°C/1 year