

Antibodies Product Datasheet

Product Name:	FAPP2 Antibody
Catalog Number:	C48718
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FAPP2
Concentration:	1 mg/ml
Gene Name:	PLEKHA8 FAPP2
Protein Name:	Pleckstrin homology domain-containing family A member 8 (PH domain-containing family A member 8) (Phosphatidylinositol-four-phosphate adapter protein 2) (FAPP-2) (Phosphoinositol 4-phosphate adapter protein 2) (hFAPP2) (Serologically defined breast cancer
Other Name:	
Human Swiss Prot No.:	Q96JA3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FAPP2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year