

Antibodies Product Datasheet

Product Name:	Cystatin B Antibody
Catalog Number:	C0330
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Stefin B. AA range:49-98
Concentration:	1 mg/ml
Gene Name:	CSTB
Protein Name:	Cystatin-B
Other Name:	CSTB; CST6; STFB; Cystatin-B; CPI-B; Liver thiol proteinase inhibitor; Stefin-B
Human Swiss Prot No.:	P04080
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cystatin B Polyclonal Antibody detects endogenous levels of Cystatin B protein.
Storage and Stability:	-20°C/1 year