

## Antibodies Product Datasheet

<b>Product Name:</b>	HNF-4 $\alpha$ / $\gamma$ Antibody
<b>Catalog Number:</b>	C10596
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human HNF4 alpha/gamma. AA range:91-140
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	HNF4A/HNF4G
<b>Protein Name:</b>	Hepatocyte nuclear factor 4-alpha/gamma
<b>Other Name:</b>	HNF4A; HNF4; NR2A1; TCF14; Hepatocyte nuclear factor 4-alpha; HNF-4-alpha; Nuclear receptor subfamily 2 group A member 1; Transcription factor 14; TCF-14; Transcription factor HNF-4; HNF4G; NR2A2; Hepatocyte nuclear factor 4-gamma; HNF-4-ga
<b>Human Swiss Prot No.:</b>	P41235/Q14541
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	HNF-4 $\alpha$ / $\gamma$ Polyclonal Antibody detects endogenous levels of HNF-4 $\alpha$ / $\gamma$ protein.
<b>Storage and Stability:</b>	-20°C/1 year