

Antibodies Product Datasheet

Product Name:	TPH1 (Phospho-Ser260) Antibody
Catalog Number:	A1012
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Tryptophan Hydroxylase around the phosphorylation site of Ser260. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	TPH1
Protein Name:	Tryptophan 5-hydroxylase 1
Other Name:	TPH1; TPH; TPRH; TRPH; Tryptophan 5-hydroxylase 1; Tryptophan 5-monooxygenase 1
Human Swiss Prot No.:	P17752
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-TPH1 (S260) Polyclonal Antibody detects endogenous levels of TPH1 protein only when phosphorylated at S260.
Storage and Stability:	-20°C/1 year