

Antibodies Product Datasheet

Product Name:	p53(Phospho-Ser366) Antibody
Catalog Number:	C40058P
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human p53(Phospho-Ser366)
Concentration:	1 mg/ml
Gene Name:	TP53 P53
Protein Name:	p53(Phospho-Ser366)
Other Name:	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
Human Swiss Prot No.:	P04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous p53(Phospho-Ser366)
Storage and Stability:	-20°C/1 year