

Antibodies Product Datasheet

Product Name:	α -tubulin Monoclonal Antibody(8F11), Cy3 Conjugated
Catalog Number:	N1238
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF/ICCIP
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC 1:50-300, IF 1:200.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse IgG
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	TUBA1A
Protein Name:	Tubulin alpha-1A chain
Other Name:	TUBA1A
Human Swiss Prot No.:	Q71U36/P68363
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	α -tubulin Monoclonal Antibody(8F11) Cy3 Conjugated specially designed for your Immunofluorescence analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in dark.