

Antibodies Product Datasheet

Product Name:	Cyclin E1 Antibody
Catalog Number:	B0636
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cyclin E1. AA range:43-92
Concentration:	1 mg/ml
Gene Name:	CCNE1
Protein Name:	G1/S-specific cyclin-E1
Other Name:	CCNE1; CCNE; G1/S-specific cyclin-E1
Human Swiss Prot No.:	P24864
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cyclin E1 Polyclonal Antibody detects endogenous levels of Cyclin E1 protein.
Storage and Stability:	-20°C/1 year