

Antibodies Product Datasheet

Product Name:	MMP-16 Antibody
Catalog Number:	C0268
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MMP-16. AA range:551-600
Concentration:	1 mg/ml
Gene Name:	MMP16
Protein Name:	Matrix metalloproteinase-16
Other Name:	MMP16; MMPX2; Matrix metalloproteinase-16; MMP-16; MMP-X2; Membrane-type matrix metalloproteinase 3; MT-MMP 3; MTMMP3; Membrane-type-3 matrix metalloproteinase; MT3-MMP; MT3MMP
Human Swiss Prot No.:	P51512
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MMP-16 Polyclonal Antibody detects endogenous levels of MMP-16 protein.
Storage and Stability:	-20°C/1 year