

Antibodies Product Datasheet

Product Name:	RSK3 (Phospho-Thr353) Antibody
Catalog Number:	A9465
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human RSK3 (Phospho Thr353)
Concentration:	1 mg/ml
Gene Name:	RPS6KA2 MAPKAPK1C RSK3
Protein Name:	RSK3 (Phospho Thr353)
Other Name:	Ribosomal protein S6 kinase alpha-2 (S6K-alpha-2) (EC 2.7.11.1) (90 kDa ribosomal protein S6 kinase 2) (p90-RSK 2) (p90RSK2) (MAP kinase-activated protein kinase 1c) (MAPK-activated protein kinase 1c) (MAPKAP kinase 1c) (MAPKAPK-1c) (Ribosomal S6 kinase 3
Human Swiss Prot No.:	Q15349
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse RSK3 (Phospho Thr353)
Storage and Stability:	-20°C/1 year