

Antibodies Product Datasheet

Product Name:	POLE3 Antibody
Catalog Number:	C48850
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human POLE3
Concentration:	1 mg/ml
Gene Name:	POLE3 CHRAC17
Protein Name:	DNA polymerase epsilon subunit 3 (EC 2.7.7.7) (Arsenic-transactivated protein) (AsTP) (Chromatin accessibility complex 17 kDa protein) (CHRAC-17) (HuCHRAC17) (DNA polymerase II subunit 3) (DNA polymerase epsilon subunit p17)
Other Name:	
Human Swiss Prot No.:	Q9NRF9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of POLE3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year